<b>UM</b>	<b>JABATAN KIMIA</b>	<b>Kategori</b>	<b>Arahan Kerja</b>
	<b>PROSEDUR PENGENDALIAN PERALATAN FT-NMR ECX 400MHz SYSTEM JOEL BAGI ANALISIS SAMPEL</b>	<b>Dokumen No.</b>	<b>AK/JK.P/C/02</b>
		<b>Tarikh   Semakan</b>	<b>27.07.2018   1</b>
		<b>Mukasurat:</b>	<b>1 daripada 10</b>

## 1. TUJUAN

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
Delta console window and file browser window

<b>UM</b>	<b>JABATAN KIMIA</b>	<b>Kategori</b>	<b>Arahan Kerja</b>
	<b>PROSEDUR PENGENDALIAN PERALATAN FT-NMR ECX 400MHz SYSTEM JOEL BAGI ANALISIS SAMPEL</b>	<b>Dokumen No.</b>	<b>AK/JK.P/C/02</b>
		<b>Tarikh   Semakan</b>	<b>27.07.2018   1</b>
		<b>Mukasurat:</b>	<b>10 daripada 10</b>

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#### 4. PROSEDUR BERKAITAN

Instrument Usage/Service Application Procedure GP-I-001

#### 5. RUJUKAN SEMAKAN

SEMAKAN	DISEDIAKAN	DISEMAK	TARIKH	CATATAN
1	Dara Fiona Mohamad  Sugakumar A/L Varuthan	Fateh Ngaliman	27.07.2018	Isu kali pertama

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>1 of 10</b>

## 1. PURPOSE

The purpose of this procedure established is to provide clear instructions for steps or orders on how to run sample analysis using JEOL NMR ECX 400MHz SYSTEM JOEL equipment.


## 2. SCOPE

Applied for user from Chemistry Department.

## 3. INSTRUCTIONS PROCEDURE


### 3.1 Delta STARTUP AND CONNECTION

- a) Air Compressor ON mode.
- b) Type Ctrl + Alt + Delete at a time.  
Username: delta  
Password: delta (displayed as asterisks (\*\*\*\*\*) ).
- c) Double-click the Delta icon  on the desktop.
- d) The Delta console window opens.

**Delta console window**


- e) Click the circled button  indicated in the figure below.
- f) The “Spectrometer control” window opens.


**SPECTROMETER CONTROL WINDOW**

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>2 of 10</b>

g) Select a spectrometer displayed in the “Spectrometer Control ” window.

h) Clicking the connection button  switches to the “User Authentication” window.

**User authentication window**


i) Type the user name and the password.

- Name delta
- Password: delta ( displayed as asterisks ( \*\*\*\*\* )

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j) Click the  button.

k) The window view changes as shown below.


**After completion of user authentication**

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	Category	Work Instruction
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	Document No.	AK/JK.P/C/01-BI
		Date   Revision	31.07.2018   1
		No. of pages:	3 of 10


### 3.2 SAMPLE PREPARATION

- a) Mount the sample tube on the sample rotor and the holder.

**Sample mounting**

- b) Set the sample tube set in the sample rotor and holder so that it floats in the SCM.  
#When using the auto-sample changer, place the tube in the slot.

**Setting the sample tube unit on the SCM**

***IMPORTANT NOTICE :***

- Before making the sample float in the SCM, check that floating air emerges.
- Check that no sample is inside the SCM.
- Do not insert a bare glass tube or the empty rotor and holder.

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>4 of 10</b>


### 3.3 PREPARATION FOR MEASUREMENT

#### (a) Creating Sample Definitions

A sample definition is:

Conditions of a sample to be measured.

1. Click the  button to create a new sample definition ( see the figure below ).

Creating sample definitions

Details of the parameters displayed in the sample definitions:

No.	Sample Name	Solvent	Slot	Kind	Shared	Verified	Error	Owner	Last Load
1		NONE	0	Liquids				delta	Never

Item	Description
	Add a sample definition
	Delete a sample definiton
<b>Sample name</b>	Arbitrary ( used as a part of a saved name )
<b>Solvent</b>	Select a solvent. Select a sample slot ( if an auto sample changer is integrated ).
<b>Slot</b>	Select a sample slot ( if an auto sample changer is integrated ).


<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>5 of 10</b>

(b) Creating a Job

**A job is:**


**A process of configuring conditions for performing an NMR measurement.**

1. Click the  button displayed in the Spectrometer Control Window.

**SPECTROMETER CONTROL WINDOW**

2. The window automatically switches to the Job tab.


**Jobs tab**

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>6 of 10</b>


### 3.4 1H MEASUREMENT

- a) Click the  button to select a pulse sequence.

Selecting a pulse sequence

- b) Select basic / proton.jxp. ( If not found, search by double clicking  ).

- c) Select the pulse sequence to use, and then click the  button.
- d) The measurement is started by clicking the  button.

**Note:**

The sample preparation, sample definition creation, job creation and measurement have been

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>7 of 10</b>

described.

These are basic procedures for performing measurements using Delta V5.


This procedure should be practiced in order to achieve a smooth measurement process.

### 3.5 <sup>1</sup>H DATA PROCESSING

When the NMR measurement is completed, the measurement result is displayed on the screen.


#### a) Peak Picking

1. Select the  button in the toolbar to determine the threshold.

#### i) Automatic detection

Peak picking is performed by clicking the  button.

After peak picking

## 2. Integration

### i. Automatic detection

1. An integral curve can be automatically detected by clicking the  button.


### ii. Adding and deleting a selected integral curve

1. To draw an integral curve separately, select the  button in the toolbar.

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>8 of 10</b>

2. An integral curve can be created by dragging a selected part on the X axis.


**Dragging right and left**

### iii. Saving processed data

1. Select **File** from the pull-down in the 1D Processor window.
2. The data can be saved by using **Save as**.



### 3.6 PRINTING DATA

a) Clicking the  button opens the printing option window.

### 3.7 <sup>13</sup>C MEASUREMENT

#### a) Measurement

- i. Create a new job or add parameters to an existing job.
  - a. Click the  button to select basic / carbon.jxp.
  - b. ( If not found, search by double clicking  ).
- ii. Change the parameters as needed.
- iii. The measurements is started by clicking the  button.

#### b) Data processing

Perform the processing for the acquired spectrum.


**Follow the same procedure as <sup>1</sup>H.**

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>9 of 10</b>

- **Phase Correction**
- **Reference setting**
- **Peak picking**

### 3.8 DEPT MEASUREMENT

#### a) Measurement

- i. Create a new job or add parameters to an existing job.
  - a. Click the  button to select basic / dept.jxp.
  - b. ( If not found, search by double clicking  ).
- ii. Change the parameters as needed.
 


[Pulse] tab: selection\_angle=135[deg]

The parameter above can be changed to 45[deg], 90[deg], or 135[deg].
- iii. The measurement is started by clicking the  button.

### 3.9 DATA SAVE LOCATION

#### a) Opening Past Data

- i. Click the  button in the [delta] console window.

**Delta console window and file browser window**

<b>UM</b>	<b>CHEMISTRY DEPARTMENT</b>	<b>Category</b>	<b>Work Instruction</b>
	<b>PROSEDURE FOR EQUIPMENT HANDLING FT-NMR ECX 400MHz SYSTEM JOEL FOR SAMPLE ANALYSIS</b>	<b>Document No.</b>	<b>AK/JK.P/C/01-BI</b>
		<b>Date   Revision</b>	<b>31.07.2018   1</b>
		<b>No. of pages:</b>	<b>10 of 10</b>

- ii. The File Browser window appears.
- iii. Click the  button.
- iv. The selected data is opened by clicking the  button.  
When the saving location for a certain data is changed, specify the location separately. The data once displayed on the screen is saved in the  folder.

#### 4. RELATED PROCEDURES

Instrument Usage/Service Application Procedure GP-I-001

#### 5. REFERENCE CHECKING

REVISION	PREPARED BY	CHECKED BY	DATE	REMARK
1	Sugakumar A/L Varuthan Dara Fiona Mohamad	Fateh Ngaliman	31.07.2018	1 <sup>st</sup> issue